

# Science. Applied to Life.™


Curing guidelines for using bulk fill resin composites, according to the Halifax Consensus Statement\*

3M<sup>™</sup> Elipar<sup>™</sup> DeepCure-S LED Curing Light 3M<sup>™</sup> Elipar<sup>™</sup> DeepCure-L LED Curing Light









### Before you start ...


- Measure the depth of the preparation
- Do not exceed product maximum increment thickness
- Follow manufacturer's recommended guidelines for light curing
- Incremental filling and light curing each increment separately may be required for large restorations
- Refer to the Instructions for Use for more details

# Before you fill a Class II ...

 Ensure you have a good matrix contour and adaption, especially at the bottom of the proximal box

# When you fill ...

- Fill carefully to avoid trapping air
- Some bulk fill composite resins should be covered with a conventional composite to improve esthetics, reduce wear and to help create occlusal anatomy



# 5 A B C



#### When you light cure ...

- Short exposure times using high-output curing lights may compromise the properties of some resins
- Cure each increment of composite according to product manufacturer's guidelines
  - A) Keep light tip close to the restoration
  - B) Avoid creating shadows
  - C) Use multiple exposures to fully cover the restoration
- Light cure from the buccal and lingual directions after removing the matrix
- Avoid overheating the tooth and gingiva (Directing high-volume suction across tooth during light curing helps)

Steps 1 through 6 are general guidelines from the Halifax Consensus Statement\* and may not be applicable to every curing light. For Elipar DeepCure-S and Elipar DeepCure-L LED Curing Lights, please refer to Instructions for Use. Resource courtesy of Dr. Richard Price, professor and head of fixed prosthodontics, department of dental clinical sciences, faculty of dentistry, Dalhousie University. Email: rbprice@dal.ca

\*Source: Halifax Consensus Statement from the 2014 Symposium on Light Curing in Dentistry, Dalhousie University, Halifax, Canada. First published in the Journal of Canadian Dental Association, j can dent assoc 2014;80:e61.

# Choose from two models, both with identical technical performance.

# 3M™ Elipar™ DeepCure-S LED Curing Light



# 3M™ Elipar™ DeepCure-L LED Curing Light



#### **Technical Performance Data for Both Models**

Wavelength	430–480 nm
Light intensity	1,470 mW/cm² (-10%/+20%)
Power supply	Lithium-ion battery
	Approx. 120 min. battery runtime (~720 10-sec. cures) with constant light output regardless of battery charge
Operation	Intuitive two-button and single-mode operation
	Pre-set cure times: 5, 10, 15 and 20 seconds, continuous mode (120 sec.) and tack-cure mode
Curing time	Refer to material instructions; 10 sec. for many composites
Light guide	10 mm, black coated, autoclavable, optimal intraoral reach due to user- and patient-friendly geometry

Ordering Information—Stainless Steel		
Item #	Product Information	
76976	3M™ Elipar™ DeepCure-S LED Curing Light Contains: Handpiece (Cordless); Charging Base; Li-ion Battery; 10mm Light Guide; Eye Shield	
76981	3M™ Elipar™ DeepCure-S Light Guide, 10mm	
76984	3M™ Elipar™ DeepCure Eye Shield	
76985	3M™ Elipar™ DeepCure-S Rechargeable Li-ion Battery	

Ordering Information—Lightweight		
Item #	Product Information	
76973	3M™ Elipar™ DeepCure-L LED Curing Light Contains: Handpiece (built-in Li-ion battery); Universal Power Supply with 5 adaptors; 10mm Light Guide; Eye Shield; 3 Curing Discs	
76983	3M™ Elipar™ DeepCure-L Light Guide, 10mm	
76984	3M™ Elipar™ DeepCure Eye Shield	
76965	3M™ Elipar™ DeepCure-L Curing Discs (5 pcs.)	

# www.3M.com



3M Oral Care 2510 Conway Avenue St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249 Web 3M.com/dental 3M Deutschland GmbH Location Seefeld 3M ESPE • ESPE Platz 82229 Seefeld • Germany Info3MESPE@mmm.com www.3MESPE.com

3M, Elipar and ESPE are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2018. All rights reserved.