

# Treatment of a deep carious lesion of second lower premolar with 3M™ Filtek™ One Bulk Fill Restorative

Clinical dentistry and photography by: Dr. Filippo Dini, Viareggio, Italy

## About the Case:

Male patient 47 years old.

Patient reported pain when drinking cold beverages. The x-ray showed the extension of caries close to pulp. Endodontic treatment was required.

## Challenge:

The deep cervical margin and the preservation of tissue for the post-endo restoration posed challenges in this case. Excellent adaptation between the material and cavity walls is critical.

## Step-by-Step:



Fig. 1: Initial Situation: Second lower right premolar with a deep carious lesion restoration that required endodontic treatment. Image is showing the final preparation after removing the temporary filling.



Fig. 2: Placement and adaptation of the sectional matrix.



Fig. 3: Enamel is etched for 15 seconds using 3M™ Scotchbond™ Universal Etchant for a selective etch approach, followed by rinsing and drying.



Fig. 4: 3M™ Scotchbond™ Universal Adhesive is applied and scrubbed into the surface for 20 seconds.



Fig. 5: After air drying for approximately 5 seconds, the adhesive is light cured for 10 seconds with 3M™ Elipar™ DeepCure-S LED Curing Light.



Fig. 6: 3M™ Filtek™ One Bulk Fill Restorative, in shade A3, is directly placed into the cavity in a single increment and then cured according to the Instructions for Use.

# Treatment of a deep carious lesion of second lower premolar with 3M™ Filtek™ One Bulk Fill Restorative (continued)



Fig. 7: Following finishing, to achieve the final polish of the restoration, the 3M™ Sof-Lex™ Pre-Polishing Spiral is used, as a first step, on a moist surface.



Fig. 8: A final high-gloss polish is created using the 3M™ Sof-Lex™ Diamond Polishing Spiral on a moist surface.



Fig. 9: Final restoration is very natural-looking and esthetic.

## The 3M Difference:

In deep and tight cavities, incremental layering can be a timely challenge and lead to poor marginal adaptation. With its easy handling and excellent adaptation, 3M™ Filtek™ One Bulk Fill Restorative enables one-step placement up to 5 mm, without compromising the esthetics of posterior restorations.

Refer to Instructions for Use (IFU) for complete product information.

# 3M

3M Oral Care  
2510 Conway Avenue  
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249  
Web 3M.com/dental

3M Canada  
Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
1-888-363-3685

3M, Elipar, Filtek, Scotchbond and Sof-Lex are trademarks of 3M Company. Used under license in Canada. All other trademarks are owned by other companies. Please recycle. Printed in U.S.A. © 3M 2016. All rights reserved.