

Kupplung	Kompatible Airscaler			
	KaVo	W&H	NSK	VDW Empfehlung
KaVo Multiflex 	SONICflex LUX 2003 LX SONICflex 2003 SONICflex 2000 N	ZA-55 LM ZA-55 M	TI-Max S970 KL	EDDY Starter Kit V04 1443 000 000 
Sirona 		ZA-55 LS	TI-Max S970 SL	TI-Max S970 SL oder ZA-55 LS
W&H Roto-Quick 		ZA-55 L ZA-55		ZA-55 L
NSK 			TI-Max S970 TI-Max S970 L	TI-Max S970 L

nicht kompatible Airscaler

KaVo SONICflex 2008
mit Quick-Kupplung

KAVO SONICflex quick Set 2008 LS

ohne Kupplung - Direktanschluss	KaVo	W&H	NSK	VDW Empfehlung
Ritter-Midwest 4-Loch 		ZE-55 RM	AS2000 M4	
Borden 2(3)-Loch 		ZE-55 BC	AS2000 B2/B3	