

Indications	Impression technique	Single-unit crown	Small implant (≤ 2 units)	Inlay/onlay (≤ 2 units)	Small bridge (≤ 3 units)	Veneers	Bridge (> 3 units)	Implant (> 2 units)	Pick-up impression	Functional impression
Super Quick Setting Materials										
3M™ Impregum™ Penta™ Super Quick (Medium Body)	Monophase									
3M™ Impregum™ Penta™ Super Quick (Medium or Heavy Body) 3M™ Impregum™ Super Quick (Light Body)	1-step	✓	✓	✓	✓					
Quick Setting Materials										
3M™ Impregum™ Penta™ Soft Quick (Medium Body)	Monophase									
3M™ Impregum™ Penta™ H DuoSoft™ Quick (Heavy Body) 3M™ Impregum™ L DuoSoft™ Quick (Light Body)	1-step	✓	✓	✓	✓	✓				
Regular Setting Materials										
3M™ Impregum™ Penta™ Soft (Medium Body)	Monophase									
3M™ Impregum™ Penta™ H DuoSoft™ (Heavy Body) 3M™ Impregum™ Penta™ L DuoSoft™ (Light Body) 3M™ Impregum™ Garant™ L DuoSoft™ (Light Body)	1-step					✓	✓	✓	✓	✓
3M™ Impregum™ Penta™ (Medium Body)	Monophase					✓	✓	✓	✓	✓