



3M Temporization Procedure	Benefits for Your Procedure	Benefits for Productivity	Benefits for Patients
Mix and Dispense 3M [™] Pentamix [™] Lite Automatic Mixing Unit	 Fewer cartridge changes Less waste than hand mix Homogenous and void-free matrices 	 Fast extrusion maximizes work time, reduces stress 	 Reduces annoying waiting times More time for clinician to care for patient
Impress 3M [™] Express [™] XT VPS Impression Materials	 Smooth matrix generates the lowest smear layer No dust 	 Real putty consistency with automix delivery More speed and improved, reproducible outcomes 	 Pleasant temperature for more patient comfort Fast-setting, thermoactive material
Fabricate Temporary 3M™ Protemp™ 4 Temporization Material	 Record-breaking toughness over time and across indications Outstanding esthetics—brilliant gloss without polishing 	 Fast and easy automix procedure Fewer remakes due to fracture Fewer unplanned patient visits 	 Fewer failures—higher patient satisfaction Natural gloss that's easy to clean and comfortable to wear
3M [™] Sof-Lex [™] Finishing and Polishing System	• Flexible discs conform to the surface of the restoration	 Efficient gross reduction and contouring 	 Adjustments are fast and easy
Cement 3M [™] RelyX [™] Temp NE/E Temporary Cement	 Good bond strength Easy excess removal Easy removal of provisionals 	 Fast and easy provisional cementation 	 Low film thickness helps provide an optimal fit
Bottom Line:	Greater convenience More control	A more productive, profitable practice	Happier, more loyal patients



Ordering Information

3M[™] Protemp[™] 4 Temporization Material

Item No.	Product Informa	tion		
46953	 3M[™] Protemp[™] 4 Temporization Material Introductory Kit—A2 Includes: 1 – 50 ml cartridge, Shade A2; 16 – 3M[™] Garant[™] Mixing Tips, Blue; 1 – 3M[™] Garant[™] Dispenser 10:1 3M[™] Protemp[™] Temporization Material Refills Each includes 1 – 50ml cartridge; 16 – 3M[™] Garant[™] Mixing Tips, Blue 			
46959	Bleach	46957	A3	
46954	A1	46972	A3.5	
46956	A2	46960	B3	
71453	3M™ Garant™ M Includes: 50 Mix	lixing Tips Refill—Bl i ing Tips	ue	

3M[™] RelyX[™] Temporary Cements

Item No.	Product Information
56660	3M™ RelyX™ Temp NE Cement (non-eugenol) Includes: 1 – 30g 3M™ RelyX™ Temp NE base paste; 1 – 14g RelyX Temp NE catalyst paste; 1 mixing pad
35014	3M™ RelyX™ Temp E Cement (eugenol) Includes: 1 – 25g 3M™ RelyX™ Temp E base paste; 1 – 18g RelyX Temp E catalyst paste; 1 mixing pad

3M[™] Protemp[™] Crown Temporization Material

These pre-formed crowns are ideal for missing tooth structure. No impression is needed, and it's malleable and light curable. It comes in nine sizes and one universal shade.





3M Oral Care 2510 Conway Avenue St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249 Web 3M.com/dental 3M Deutschland GmbH Location Seefeld 3M ESPE • ESPE Platz 82229 Seefeld • Germany Info3MESPE@mmm.com www.3MESPE.com

3M[™] Express[™] XT Putty Soft VPS Impression Material

Item No. Product Information

36996 3M[™] Express[™] XT Putty Soft Intro Kit

Includes: 1 – 250 ml jar base paste; 1 – 250 ml jar catalyst paste; 2 putty spoons; 1 – 50 ml cartridge 3M™ Express™ XT Light Body; 1 – 50 ml cartridge 3M™ Express™ XT Regular Body; 10 – 3M™ Garant™ Mixing Tips, Yellow; 10 – 3M™ Garant™ Intra-oral Tips, Yellow; 1 VPS tray adhesive

36972 3M™ Express™ XT Putty Soft Refill

Includes: 1 – 250 ml jar base paste; 1 – 250 ml jar catalyst paste; 2 putty spoons

3M[™] Express[™] XT Penta[™] Putty VPS Impression Material

Item No. Product Information

36892	3M™ Express™ XT Penta™ Putty Single Pack
	Includes: 1 – 300 ml base paste; 1 – 60 ml catalyst paste

3M[™] Pentamix[™] Lite Automatic Mixing Unit and Accessories

Product Information
3M™ Pentamix™ Lite Automatic Mixing Unit
Includes: 1 – 3M™ Penta™ Cartridge for 3M™ Pentamix™ Lite
3M™ Penta™ Cartridge for 3M™ Pentamix™ Lite
3M™ Penta™ Mixing Tips—Red (50 mixing tips)

www.3M.com

3M, ESPE, Express, Garant, Penta, Pentamix, Protemp, RelyX and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2017. All rights reserved.