

# Now our most popular flowable is easier to inject.

3M<sup>™</sup> Filtek<sup>™</sup> Supreme Flowable Restorative is available in an ergonomic syringe designed for improved comfort and easier injection.\*



Virtually no bubbles or material "run-on" during dispensing gives you better control.\*

#### No more guessing

The blue syringe barrel indicates the remaining material volume. When your syringe is empty, only the white plunger shows.





Bendable cannula Easy-to-bend without kinking for better access.\* Easy to hold and inject

Triangle finger plate and plunger make it quick and easy to personalize your grip.\*

# **Key benefits**

- Virtually no bubbles or material "run-on"\*
- Easy-to-bend cannula for improved access\*
- Ergonomic design is easy to hold and inject\*
- Excellent adaptation, polish retention and wear resistance\*
- Versatile for a wide range of indications and techniques

#### Can be warmed for injection moulding

Warm up to 70°C for up to one hour, up to 25 times.\*



# **Clinical cases**

3M<sup>™</sup> Filtek<sup>™</sup> Supreme Flowable Restorative offers a smooth transition from material to tooth for aesthetic restorations.



Before

**Flowable placement** 



After

3M<sup>™</sup> Filtek<sup>™</sup> Supreme Flowable Restorative adapts easily and seals

Photos courtesy of Dr. Gunnar Reich, Munich, Germany

Photos courtesy of Dr. Giuseppe Chiodera, Brescia, Italy

the margin.




3M.co.uk/Dental

Before

**Flowable placement** 

After





#### References

1. Including small, non stress-bearing occlusal restorations. 2. Repair of acrylic temporary and composite resin.

\*3M Internal data

**3M Oral Care 3M United Kingdom PLC** Charnwood Campus 10 Bakewell Road Loughborough Leicestershire, LE11 5RB 0800 626 578

3M Ireland Ltd The Iveagh Building Carrickmines Park Carrickmines Dublin 18 Ireland +353 1 280 3555

Indications (Incremental up to 2 mm) Class III Class V **Minimally Invasive** Dentistry (MID)<sup>1</sup> Base Liner **Repair of Indirect** (small defects) Repair<sup>2</sup> Undercut Pit & Fissure Block-out Sealant

Refer to Instructions for Use for full list of indications.

### **Ordering information**



3M<sup>™</sup> Filtek<sup>™</sup> Supreme Flowable **Restorative Refills** Syringe refills include: 2 - 2g syringes; 20 dispensing tips.

SHADE	ITEM NO.
A1	6032A1
A2	6032A2
A3	6032A3
A3.5	6032A3.5
A4	6032A4
B1	6032B1
B2	6032B2
C2	6032C2
D2	6032D2
White	6032W
Extra White	6032XW
OA3	6032OA3

3M and Filtek are trademarks of the 3M company.	© 3M 2020. All
rights reserved. 13432284-20200402-0112A	